


Epidiascope Vh 2

500/1000 WATTS

Leitz
WETZLAR



**A Modern Visual-Aid for Lecture Hall
and Classroom**

330-3/Engl.

Large 6x7 $\frac{1}{4}$ " object area
Sliding object stage for projecting
sections of originals
Increased screen brilliance
New tangential cooling system
Projection distances up to 26 ft.

Epidiascope Vh 2



LEITZ episcopes and epidiascopes have proved their worth through decades of practical classroom and lecture-hall experience.

Their unusually robust construction guarantees years of faithful service with excellent optical performance.

The model Vh 2 is the successor to the well-known Vh type. Salient features of this new design are increased screen brilliance, extended object areas, and the possibility of projecting larger originals by sliding displacement of the object stage. This eliminates the need to disturb the original itself.

Leitz
WETZLAR


Episcopy

Definition and screen brilliance are dependent upon the technical arrangement of illumination and optical components in an epidiascope. The projection of episcopic (opaque) originals makes especially high demands of such an apparatus. These demands are ideally met by the illumination arrangement in the LEITZ Vh2 epidiascope, in conjunction with highly corrected LEITZ projection anastigmats.

The light-beam from the projection lamp is concentrated on the object area by means of a reflector and three illuminating mirrors so that the entire area is evenly illuminated, from corner to corner. High-speed LEITZ projection anastigmats yield clear, brilliant screen images, free from distortion and with accurate color rendition.

A quietly running tangential blower, in conjunction with a special heat filter, provides optimum cooling of the episcopic object area. Projection lamp and blower are switched on and off by means of a common switch.

The episcopic object area measures 6" (16cm) high by 7 1/4" (19cm) wide. The sturdy, horizontally adjustable object stage permits the sectional projection of originals measuring up to 7 1/4" x 12 1/4" (19 x 31cm). The objects can thus be "driven" continuously without having to be turned round and repositioned at intervals. Since the object stage is adjusted horizontally together with the original, the projection lens remains at a fixed distance from the screen. No variations in image scale or definition can therefore occur on displacing a section of the original. Moreover, the object stage can be secured vertically in any position by means of a cam lever, thus enabling originals of various thicknesses to be projected. Printed, handwritten or drawn originals, books, tables, maps or atlases, to mention but a few examples, can therefore be projected in rapid succession.



6344-32

Episcopic object area 7 1/4" x 12 1/4" (19 x 31cm) provided by sliding object stage.


Diascopy

The slide attachment is detachable, thus enabling a model originally supplied solely for episcopic projection to be converted into an epidiascope. The large stage is designed to accept slides 3 1/4" x 3 1/4" (8.5 x 8.5cm), 3 1/4" x 4" (8.5 x 10cm) and 3 1/2" x 4 3/4" (9 x 12cm). The optical illumination arrangement consists of a triple-lens condensor system.

The transition from episcopic to diascope projection is effected simply by turning a control knob.

Projection lamps

500 or 1000 watt projection lamps may be used, the latter in conjunction with a supplementary heat filter for the diascope arrangement.



Episcopic ray path



Diascopic ray path

Projection Distance and Screen Size

The focal length of the lens for episcopy and diascopy should be selected according to the space available (or desired projection distance) and screen size. Lenses are so matched in the basic epidiascopic outfit, that the episcopic and diascope images will be of comparable sizes.

Dimensions:

Base (width x length)	10 ⁵ / ₈ x 35 ⁵ / ₈ "
Height	43"
Length with diascopic lens	1 : 4 / 200mm
Length with diascopic lens	23 ¹ / ₂ "
Length with diascopic lens	1 : 4 / 250mm
Weight	45 ³ / ₈ "
	84 lbs

Outfits

		Episcopy Cat. No.	Epidiascopy Cat. No.
Episcopy Vh 2 less lamp			
with 400mm EPIS f/3.6	33 700		
with 400mm EPIS f/3.6 and 200mm DIMAR f/4			33 705
with 400mm EPIS f/3.6 and 250mm DIMAR f/4			33 715
with 500mm EPIS f/4.3	33 720		
with 500mm EPIS f/4.3 and 250mm DIMAR f/4			33 725

Design subject to alteration without notice.

Type of Projection	Projection distance in ft. (screen to face of apparatus base)				
Focal Length of Lenses	13	16	20	23	26
Episcopy	Screen size for episcopic side length 7 ¹ / ₄ " (19cm)				
400mm object area 6 ⁵ / ₁₆ x 7 ¹ / ₂ " (16 x 19cm)	6'	7' 7"	9' 2"	10' 8"	
500mm object area 6 ⁵ / ₁₆ x 7 ¹ / ₂ " (16 x 19cm)	4' 9"	6'	7' 3"	8' 4"	9' 8"
Diascopy					
Side length of slide 3 ³ / ₈ " (8.5cm)	4' 7"	6'	7'	8' 2"	
200mm Side length of slide 3 ¹ / ₂ " (10cm)	5' 7"	7'	8' 4"	9' 10"	
Side length of slide 4 ³ / ₄ " (12cm)	6' 9"	8' 6"	10' 4"	12' 2"	
Side length of slide 3 ³ / ₈ " (8.5cm)	3' 7"	4' 7"	5' 7"	6' 7"	7' 7"
250mm Side length of slide 3 ¹ / ₂ " (10cm)	4' 3"	5' 5"	6' 7"	7' 9"	9'
Side length of slide 4 ³ / ₄ " (12cm)	5' 3"	6' 7"	8'	9' 6"	11'

Accessories

	Cat. No.
Projection lamp, 500-watt (pin base); please state voltage when ordering	37 707
Projection lamp, 1000-watt / 110 volts	37 706
Projection lamp, 1000-watt / 220 volts *	37 713
Heat filter for diascopic outfit in conjunction with 1000-watt lamp	37 840
Variable transformer, 500-watt, 220 / 110 volts	37 500
Fixed transformer, 1000-watt, 220 / 110 volts	37 501
Variable transformer, 1000-watt, 220 / 110 volts	37 502
Protective dust-cover on request	

* 110 volt lamps yield more intense light than 220 volt lamps.

Leitz
WETZLAR

ERNST LEITZ GMBH WETZLAR GERMANY

Subsidiary:
Ernst Leitz (Canada) Ltd., Midland, Ontario

List **330 - 3 / Engl.**

Printed in Germany

VIII/67/FY/B