

MONOCULAR and BINOCULAR

Medical

Microscope SM


INTRODUCTION

This catalog describes our microscope Model SM. It constitutes a radical departure from conventional design; and in the short time since it appeared on the market, it has been acclaimed as the outstanding instrument for studies by the leading medical schools throughout the United States.

The specifications which are listed hereunder fill even the most exacting requirements of the microscope committee at your University.

The selection of a microscope for medical studies should be made very carefully. It is a life-time investment, and your choice of microscope must consider the benefits which superior optical as well as mechanical performance will provide in your studies. LEITZ Microscopes have been manufactured for over 120 years and have an enviable record of dependability, optical quality and forward design.

The latest technological advances in the scientific instrument field are incorporated in the design of the SM. Coarse and fine focusing motions are combined in a single operating control.




This makes focusing easier and entirely free from lost motion. The mechanism operates smoothly under all climatic conditions and independent of lubricants.


A great deal of emphasis has been placed upon the ready interchangeability of the essential parts, such as tubes, condensers, mirror and fork for external light sources, and the built-on microscope lamps.

A new improved binocular tube doubles the light transmission and insures a brillant image even with the standard substage illuminator. All medical student microscopes are equipped with widefield eyepieces.

The objective carrier is a dust-proof nosepiece operating on ball bearings. All objectives are par-focal and the high power lenses are equipped with spring loaded mounts to afford complete specimen protection as well as to prevent damage to the front lens of the objective.

The LEITZ SM Microscope is built for years of service and satisfaction. It is compact for carrying by medical students, and desirable for prolonged use in student laboratories.












TUBE "O" #512 064 \$25.00



TUBE "P" #512 065 \$46.00



TUBE "S" #512 066 \$213.00



PLAIN SQUARE STAGE #19 #512 028 \$14.00 GRADUATED MECHANICAL STAGE #26



#512 123 \$79.00

COAXIAL,



GRADUATED MECHANICAL STAGE #47 #051 322 \$94.00

COAXIAL, GRADUATED MECHANICAL STAGE #4 #512 025 \$168.00



SINGLE LENS CONDENSER #65 #512 020 \$22.00



TWO LENS CONDENSER #66 #512 021 \$36.00



TWO LENS CONDENSER #70 #051 342 \$15.00



SWING OUT CONDENSER #60: #512 085 \$51.00



PHASE CONTRAST CONDENSER #402A 513 140 \$199.00



IMMERSION DARKFIELD CONDENSER #82 #513 109 \$91.00



DRY DARKFIELD CONDENSER #84 #513 111 \$120.00




MIRROR IN MOUNT #512 078 \$12.00



MICRO DIA LAMP #514 096 \$23.00



LOW VOLTAGE LAMP WITH REGULATING TRANSFORMER #051 475 \$56.00




SM-M 0.4.10.12 P 19/65 Optics: Alla Mono




SM-M 0.4.10.13 P 19/65 Optics: Alla Mono

		LEITZ MONOCULAR MEDICAL MICROSCOPE SM-M CONSISTING OF:		
051	230	Modern stand, SM-M, made of non-corroding alloy with single knob, combined coarse and fine adjustment. Fixed cylindrical sleeve10, permitting the interchange of sleeve mounted condensers; including protective dust cover	\$	92.00
512	065	Interchangeable, inclined monocular tube "P"		46.00
512	028	Plain stage #19, with one pair of slide clips, 130 x 125mm		14.00
512	020	Single lens condenser #65, A. 0.65, with iris diaphragm, filter holder with blue and ground glass filters, in cylindrical sleeve mount		22.00
512	043	Quadruple revolving objective nosepiece 0.4, on ball bearing races		35.00
512	078	Plane and concave mirror12		12.00
512	210	LEITZ Microscope SM-M 0.4.10.12 P 19/65 as described above	\$	221,00
		ODDICAL EQUIDADAM ALL MONO		
E10	004	OPTICAL EQUIPMENT Alla MONO	ф	25 00
	004	Achromatic dry objective, 10x/0.25, focal length 16.3mm, free working distance 5.7mm	\$	35.00
519	003	Achromatic dry objective, $40x/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount		54.00
519	221	Huygens eyepiece, single H10x, field of view 14mm		12.00
510	047	LEITZ Monocular Medical Microscope SM-M 0.4.10.12 P $19/65$ + Alla Mono, complete as described above	\$	322.00
510	048	LEITZ Monocular Medical Microscope SM-M 0.4.10.13 P 19/65 + A11a Mono, complete as described above, however, equipped with Micro-dia Lamp 15 watts, 110 volts, in place of the mirror	\$	333.00
		OPTIONAL EQUIPMENT		
519	036	Achromatic dry objective, 3.5x/0.10, focal length 31.6mm, free working distance 23mm	\$	32.00
519	006	Achromatic dry objective, $25 \text{x}/0.50$, focal length 7.1mm, free working distance 0.92mm, with spring loaded mount		68.00
512	045	Leatherette carrying case, upright	\$	34.00
051	321	Attachable mechanical stage SK1, traversing area 75 x 26mm	\$	39.00
513	155	Attachable mechanical stage #21, with low set coaxial control knobs, traversing area $75 \times 26 \mathrm{mm}$		67.00




SM-M 0.4.11.12 P 19/601 Optics: A6a Mono




SM-M 0.4.11.13 P 19/601 Optics: A6a Mono

	LEITZ MONOCULAR MEDICAL MICROSCOPE SM-M CONSISTING OF:		
051 231	Modern stand, SM-M, made of non-corroding alloy with single knob, combined coarse and fine adjustment. Fixed slider11, accepting dovetail mount condensers; including protective dust cover	\$	97.00
512 065	Interchangeable, inclined monocular tube "P"		46.00
512 028	Plain stage #19, with one pair of slide clips, 130 x 125mm		14.00
512 085	Swing-out condenser with aperture diaphragm and interchangeable aspherical top element N.A. $0.90,\ \#601$		51.00
512 043	Quadruple revolving objective nosepiece 0.4, on ball bearing races		35.00
512 078	Plane and concave mirror12		12.00
512 214	LEITZ Microscope SM-M 0.4.11.12 P 19/601 as described above	\$	255.00
	OPTICAL EQUIPMENT A6a MONO		
519 036	Achromatic dry objective, $3.5x/0.10$, focal length 31.6 mm, free working distance 23 mm	\$	32.00
519 004	Achromatic dry objective, $10x/0.25$, focal length 16.3mm, free working distance 5.7mm		35.00
519 003	Achromatic dry objective, $40x/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount		54.00
519 221	Huygens eyepiece, single H10x, field of view 14mm	_	12.00
510 051	LEITZ Monocular Medical Microscope SM-M 0.4.11.12 P $19/601$ + A6a Mono, complete as described above	\$	388.00
510 052	LEITZ Monocular Medical Microscope SM-M 0.4.11.13 P $19/601$ + A6a Mono, complete as described above, however, equipped with Micro-dia Lamp 15 watts, 110 volts, in place of the mirror	\$	399.00
	OPTIONAL EQUIPMENT		
519 060	Achromatic oil immersion objective, Oil $100x/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount	\$	92.00
512 045	Leatherette carrying case, upright	\$	34.00
051 321	Attachable mechanical stage SK 1, traversing area 75 x 26mm	\$	39.00
513 155	Attachable mechanical stage #21, with low set coaxial control knobs, traversing area $75 \times 26 \mathrm{mm}$		67.00



SM-M 0.4.5.13 P 26/70 Optics: A2n Mono

		LEITZ MONOCULAR MEDICAL MICROSCOPE SM-M CONSISTING OF:		
051	235	Modern stand, SM-M, made of non-corroding alloy with single knob, combined coarse and fine adjustment. Rack and pinion with dovetail slider5, for the adjustment in the height of condenser; including protective dust cover	\$	120.00
512	065	Interchangeable, inclined monocular tube "P"		46.00
512	123	Mechanical stage #26, with scales and verniers, traversing area 76 x 50mm (3" x 2")		79.00
051	342	Two lens condenser $\#70$, A. 1.20, with iris diaphragm, filter holder with blue and ground glass filters, in dovetail mount		15.00
512	043	Quadruple revolving objective nosepiece 0.4, on ball bearing races		35.00
514	096	Micro-dia lamp13, 15 watt, for direct connection to 110 volts A.C.	_	23.00
		LEITZ Microscope SM-M 0.4.5.13 P 26/70 as described above	\$	318.00
		OPTICAL EQUIPMENT A2n MONO		
519	036	Achromatic dry objective, $3.5x/0.10$, focal length 31.6mm, free working distance 23mm	\$	32.00
519	004	Achromatic dry objective, $10x/0.25$, focal length 16.3mm, free working distance 5.7mm		35.00
519	003	Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount		54.00
519	060	Achromatic oil immersion objective, Oil $100 x/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount		92.00
519	128	Periplanatic eyepiece, single $10x^{O}$, field of view 18mm		29.00
		LEITZ Monocular Medical Microscope SM-M 0.4.5.13 P $26/70$ + A2n Mono complete as described above	\$	560.00
		OPTIONAL		
512	045	Leatherette carrying case, upright	\$	34.00




SM-M 0.4.5.051 475 S 26/601 Optics: A2n Bino

	LEITZ BINOCULAR MEDICAL MICROSCOPE SM-M CONSISTING OF:	
051 235	Modern stand, SM-M, made of non-corroding alloy with single knob, combined coarse and fine adjustment. Rack and pinion with dovetail slider5, for the adjustment in the height of condenser; including protective dust cover	\$ 120.00
512 066	Interchangeable, inclined binocular tube "S" (tube factor 1.25x)	213.00
512 123	Mechanical stage #26, with scales and verniers, traversing area 76 x 50mm (3" x 2")	79.00
512 085	Swing-out condenser with aperture diaphragm and interchangeable aspherical top element N.A. 0.90, $\#601$	51.00
512 043	Quadruple revolving objective nosepiece 0.4, on ball bearing races	35.00
051 475	Low voltage lamp 6.3 volts, 0.8 amps, 5 watt with adjustable aspherical condenser, blue and ground glass filters and regulating transformer (mounted to the microscope by means of a bayonet lock)	56.00
	LEITZ Microscope SM-M 0.4.5.051 475 S 26/601 as described above	\$ 554.00
	OPTICAL EQUIPMENT A2n BINO	
519 036	Achromatic dry objective, 3.5x/0.10, focal length 31.6mm, free working distance 23mm	\$ 32.00
519 004	Achromatic dry objective, $10x/0.25$, focal length 16.3mm, free working distance 5.7mm	35.00
519 003	Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount	54.00
519 060	Achromatic oil immersion objective, Oil $100 \mathrm{x}/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount	92.00
519 129	Periplanatic eyepieces, paired $10x^0$, field of view 18mm	58.00
	LEITZ Binocular Medical Microscope SM-M 0.4.5.051 475 S $26/601$ + A2n Bino complete as described above	\$ 825.00
	OPTIONAL	
512 045	Leatherette carrying case, upright	\$ 34.00




SM-M 0.4.5.051 475 S 47/601 Optics: A2n Bino

		LEITZ BINOCULAR MEDICAL AND LABORATORY MICROSCOPE SM-M CONSISTING OF:		
051	235	Modern stand, SM-M, made of non-corroding alloy with single knob, combined coarse and fine adjustment. Rack and pinion with dovetail slider5, for the adjustment in the height of condenser; including protective dust cover	\$	120.00
512 (066	Irterchangeable, inclined binocular tube "S" (tube factor 1.25x)		213.00
051	322	Mechanical stage #47, 135 x 135mm, with scales and verniers, low set operating knobs, traversing area $76 \times 50 \text{mm}$		94.00
512 (085	Swing-out condenser with aperture diaphragm and interchangeable aspherical top element N.A. 0.90, $\#601$		51.00
512 (043	Quadruple revolving objective nosepiece 0.4, on ball bearing races		35.00
051 4	475	Low voltage lamp 6.3 volts, 0.8 amps, 5 watt with adjustable aspherical condenser, blue and ground glass filters and regulating transformer (mounted to the microscope by means of a bayonet lock)		56.00
		LEITZ Microscope SM-M 0.4.5.051 475 S 47/601 as described above	\$	569.00
		OPTICAL EQUIPMENT A2n BINO		
519 (036	$\frac{\text{OPTICAL EQUIPMENT A2n BINO}}{\text{Achromatic dry objective, } 3.5\text{x}/0.10, \text{ focal length } 31.6\text{mm}, \text{ free working distance } 23\text{mm}}$	\$	32.00
519 (519 (\$	32.00 35.00
	004	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm	\$	
519 (004 003	Achromatic dry objective, $3.5x/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10x/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40x/0.65$, focal length 4.46mm, free working distance 0.67mm,	\$	35.00
519 (519 (004 003 060	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working	\$	35.00 54.00
519 (519 (519 (004 003 060	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount	\$	35.00 54.00 92.00
519 (519 (519 (004 003 060	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount Periplanatic eyepieces, paired 10x^0 , field of view 18mm LEITZ Binocular Medical and Laboratory Microscope SM-M 0.4.5.051 475 S 47/601 + A2n	_	35.00 54.00 92.00 58.00
519 (519 (519 (004 003 060	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount Periplanatic eyepieces, paired 10x^0 , field of view 18mm LEITZ Binocular Medical and Laboratory Microscope SM-M 0.4.5.051 475 S 47/601 + A2n	_	35.00 54.00 92.00 58.00
519 (519 (519 (004 003 060 129	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm , free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm , free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm , free working distance 0.67mm , with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm , free working distance 0.13mm , with spring loaded mount Periplanatic eyepieces, paired 10x^0 , field of view 18mm LEITZ Binocular Medical and Laboratory Microscope SM-M $0.4.5.051$ 475 S $47/601$ + A2n Bino complete as described above	_	35.00 54.00 92.00 58.00



SM-M 0.4.5.051 475 S 48/601 Optics: A2n Bino

		LEITZ BINOCULAR STUDENT AND LABORATORY MICROSCOPE SM-M CONSISTING OF:		
051	235	Modern stand, SM-M, made of non-corroding alloy with single knob, combined coarse and fine adjustment. Rack and pinion with dovetail slider5, for the adjustment in the height of condenser; including protective dust cover	\$	120.00
512	066	Interchangeable, inclined binocular tube "S" (tube factor 1.25x)		213.00
512	025	Mechanical stage #48, $160 \times 138 \text{mm}$, with scales and verniers, low set operating knobs on one axis, traversing area $76 \times 50 \text{mm}$ (3" x 2")		168.00
512	085	Swing-out condenser with aperture diaphragm and interchangeable aspherical top element N.A. 0.90, $\#601$		51.00
512	043	Quadruple revolving objective nosepiece 0.4, on ball bearing races		35.00
051	475	Low voltage lamp 6.3 volts, 0.8 amps, 5 watt with adjustable aspherical condenser, blue and ground glass filters and regulating transformer (mounted to the microscope by means of a bayonet lock)		56.00
		LEITZ Microscope SM-M $0.4.5.051$ 475 S $48/601$ as described above	\$	643.00
		CD-1017 TOWNS IN THE 10 THE		
		OPTICAL EQUIPMENT A2n BINO		
519	036	$\frac{OPTICAL\ EQUIPMENT\ A2n\ BINO}{Achromatic\ dry\ objective,\ 3.5x/0.10,\ focal\ length\ 31.6mm,\ free\ working\ distance\ 23mm}$	\$	32.00
	036 004		\$	32.00 35.00
519		Achromatic dry objective, 3.5x/0.10, focal length 31.6mm, free working distance 23mm	\$	
519 519	004	Achromatic dry objective, $3.5x/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10x/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40x/0.65$, focal length 4.46mm, free working distance 0.67mm,	\$	35.00
519 519 519	004 003	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working	\$	35.00 54.00
519 519 519	004 003 060	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount	\$	35.00 54.00 92.00
519 519 519	004 003 060	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount Periplanatic eyepieces, paired 10x^0 , field of view 18mm LEITZ Binocular Student and Laboratory Microscope SM-M 0.4.5.051 475 S 48/601 + A2n	_	35.00 54.00 92.00 58.00
519 519 519	004 003 060	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount Periplanatic eyepieces, paired 10x^0 , field of view 18mm LEITZ Binocular Student and Laboratory Microscope SM-M 0.4.5.051 475 S $48/601 + \text{A2n}$ Bino complete as described above	_	35.00 54.00 92.00 58.00
519 519 519 519	004 003 060	Achromatic dry objective, $3.5 \text{x}/0.10$, focal length 31.6mm, free working distance 23mm Achromatic dry objective, $10 \text{x}/0.25$, focal length 16.3mm, free working distance 5.7mm Achromatic dry objective, $40 \text{x}/0.65$, focal length 4.46mm, free working distance 0.67mm, with spring loaded mount Achromatic oil immersion objective, Oil $100 \text{x}/1.30$, focal length 1.9mm, free working distance 0.13mm, with spring loaded mount Periplanatic eyepieces, paired 10x^0 , field of view 18mm LEITZ Binocular Student and Laboratory Microscope SM-M 0.4.5.051 475 S 48/601 + A2n	_	35.00 54.00 92.00 58.00



SM-M 0.4.5.051 475 S 48/402a Optics: C 10 Bino

		LEITZ BINOCULAR STUDENT & LABORATORY PHASE MICROSCOPE SM-M CONSISTING OF:		
051	235	Modern stand, SM-M, made of non-corroding alloy with single knob, combined coarse and fine adjustment. Rack and pinion with dovetail slider5, for the adjustment in the height of condenser; including protective dust cover	\$	120.00
512	066	Interchangeable, inclined binocular tube "S" (tube factor 1.25x)		213.00
512	025	Mechanical stage #48, $160 \times 138 \mathrm{mm}$, with scales and verniers, low set operating knobs on one axis, traversing area $76 \times 50 \mathrm{mm}$ (3" x 2")		168.00
513	140	Phase contrast condenser, PHACO #402a, with aperture diaphragm, swing-out upper element Achr. 0.90 and revolving disc containing three annular diaphragms, in dovetail mount		199.00
512	043	Quadruple revolving objective nosepiece 0.4, on ball bearing races		35.00
051	475	Low voltage lamp 6.3 volts, 0.8 amps, 5 watt with adjustable aspherical condenser, blue and ground glass filters and regulating transformer (mounted to the microscope by means of a bayonet lock)		56.00
		LEITZ Microscope SM-M 0.4.5.051 475 S 48/402a as described above	\$	791.00
		OPTICAL EQUIPMENT C10 BINO		
519	165	Achromatic dry phase objective, PHACO 10x/0.25, free working distance 7.7mm	\$	80.00
519	236	Achromatic dry phase contrast objective, PHACO $25x/0.50$, focal length 7.2mm, free working distance 1.06mm, with spring loaded mount		144.00
519	166	Achromatic dry phase objective, PHACO $40 \text{x}/0.65$, free working distance 0.71mm, with spring loaded mount		102.00
519	167	Achromatic oil immersion phase objective, PHACO Oil $100x/1.30$, free working distance 0.17mm, with spring loaded mount		154.00
519	142	Periplanatic widefield eyepieces, paired GF 10x, field of view 18mm		78.00
513	123	Magnifier for focusing the phase ring		34.00
		LEITZ Binocular Student and Laboratory Microscope SM-M 0.4.5.051 475 S $48/402a$ + C10 Bino 513 123, complete as described above	\$1	, 383.00
		OPTIONAL		
512	045			
012	0.10	Leatherette carrying case, upright	\$	34.00

OBJECTIVES FOR TRANSMITTED-LIGHT INVESTIGATION IN BRIGHT-DARKFIELD

Designation of Objectives	Magnification/Aperture	Type of eyepiece	Cover glass correc- tion	Focal length mm	Free work- ing dis- tance mm	Catalog #	Price
ACHROMATIC DRY OBJECTIVES	$\begin{array}{c} 2.5 \text{x}/0.07 \\ 3.2 \text{x}/0.12 \\ 3.5 \text{x}/0.10 \\ 6 \text{x}/0.18 \\ 10 \text{x}/0.25 \\ 13 \text{x}/0.40 \\ 25 \text{x}/0.50 \\ 40 \text{x}/0.65 \\ 63 \text{x}/0.85 \\ \text{Iris } 63 \text{x}/0.85 \end{array}$	P H H H P P	D O D O D O D D D D D D D D D D D D D D	56.8 39.8 31.6 23.1 16.3 13.3 7.1 4.46 2.9 2.9	14 35 23 17 5.7 3.4 0.92 0.67 0.29	519 044 519 005 519 036 519 001 519 004 519 237 519 006 519 003 519 002 519 047	Discontinued Discontinued 32.00 34.00 35.00 59.00 68.00 54.00 77.00 111.00
ACHROMATIC IMMERSION OBJECTIVES We = WATER	Oil & We 22x/0.65 We 90x/1.20 Oil 100x/1.30 Iris Oil 100x/1.30-1.10	P P P	D O D D ⁴ D	8.1 2.1 1.9 1.9	0.32 0.09 0.13 0.13	519 059 519 171 519 060 519 045	104.00 154.00 92.00 154.00
FLUORITE DRY OBJECTIVES	Fl 40x/0.85	P	D 1	4.3	0.38	519 025	151.00
FLUORITE IMMERSION OBJECTIVES	Fl Oil 54x/0.95 Fl Oil 95x/1.32 Iris Fl Oil 95x/1.32	P P P	D O D ⁴ D	3.4 2.0 2.0	0.22 0.15 0.15	519 027 519 026 519 046	160.00 203.00 281.00
PLANO ACHROMATIC OBJECTIVES (NPL) FOR FIELD OF VIEW UP TO 18MM \$\phi\$	NPL 6.3x/0.20 NPL 16x/0.40 NPL 25x/0.50 NPL 40x/0.65 NPL Oil 100x/1.30	P P P P	D O D D D D	23.6 10.8 6.99 4.53 1.73	2.0 0.50 0.38 0.15 0.21	519 245 519 246 519 247 519 248 519 249	104.00 107.00 119.00 140.00 240.00

OBJECTIVES TRANSMITTED LIGHT, PHASE CONTRAST

Designation of Objectives	Magnification/	Aperture	Type of eyepiece	Cover glass correc- tion	Contrast	Free work- ing dis- tance mm	Catalog #	Price
PHASE CONTRAST, PHACO WITH ZERNIKE CONDENSER	PHACO PHACO PHACO PHACO Oil	10x/0.25 25x/0.50 40x/0.65 100x/1.30	P P P	D O O D	Normal Normal Normal	7.7 1.06 0.71 0.17	519 165 519 236 519 166 519 167	80.00 144.00 102.00 154.00
PLANO-PHASE CONTRAST, PHACO-NPL WITH ZERNIKE CONDENSER	PHACO-NPL PHACO-NPL PHACO-NPL PHACO-NPL PHACO-NPL Oil	10x/0.25 16x/0.40 25x/0.50 40x/0.65 100x/1.30	P P P	DO DO D	Normal Normal Normal Normal	0.50 0.38 0.15	519 250 519 251 519 252 519 253 519 254	149.00 154.00 176.00 197.00

- 1) D : with cover glass D = 0.17 (coverglass thickness should be observed to within \pm 0.05mm)
 - O: without coverglass
 - DO: can be used with or without coverglass
 - D: coverglass thickness 0.17mm should be observed accurately to within ± 0.01mm, or should be accurately set with the correction mount where it
- 20 varies from this values.

- 3) H = use Huygens eyepiece, P = use PERIPLAN or widefield eyepieces
- 4) These oil immersion objectives may also be used for uncovered subjects (smear preparations without coverglass); the negligible reduction of image quality can be ignored.

All objectives from 3.5/0.10 are parfocal on the nosepiece.

EYEPIECES

HUYGENS EYEPIECES 23.2mm φ						PERIPLANATIC WIDEFIELD EYEPIECES 23.2mm ø						ım ø
Magnification	Field of view	C.	ATALOG	NUMBER		Mag	nification	Field of view	(CATALOG	NUMBE	R
	mm	Single	Price	Pair	Price			mm	Single	Price	Pair	Price
* 6x * 6x H 6.3x H 6.3xM H 10x	19 17 18 18 14	519 034 519 180 519 182 519 221 TIC EYEF	11.00 16.00 24.00 12.00	519 030 519 181 519 183 519 222	22.00 32.00 40.00 24.00	GF GF GF	10x 10xM 12.5x 12.5xM 16x 25x 25xM	18 18 18 18 15 10	519 137 519 126 519 051 519 055 519 138 519 140 519 141	39.00 50.00 49.00 56.00 46.00 50.00 58.00	519 142 519 127 519 053 519 056 519 143 519 144 519 130	78.00 89.00 98.00 105.00 92.00 100.00
Magnification	Field of view			NUMBER		-	nification	Field of view			NUMBE	
	mm	Single	Price	Pair	Price			mm	Single	Price	Pair	Price
6.3x 6.3xM 8x 10x 10xM * 12x 12.5xMM 15x	18 18 16 18 14.6 13 12.5	519 185 519 188 519 139 519 128 519 170 519 219 569 059 051 982	33.00 41.00 24.00 29.00 44.00 23.00 67.00 23.00	519 186 519 187 519 028 519 129 519 035 051 980 569 058 051 983	66.00 74.00 48.00 58.00 88.00 46.00 134.00 46.00	8: 10: 10:		18 15 15	519 147 519 135 519 134	46.00 48.00 62.00	519 148 519 136 519 125	92.00 96.00 110.00

^{*} Discontinued, limited supply still available

INDIVIDUAL COMPONENTS

		TUBES	
5 1	2 064	Interchangeable, straight monocular tube "O"	\$ 25.00
51	2 065	Interchangeable, inclined monocular tube "P"	46.00
51	2 066	Interchangeable, inclined binocular tube "S" (tube factor 1.25x)	213.00
51	.3 264	Discussion tube for the simultaneous observation by two people of the microscope image with built-in mechanical arrow pointer and two bayonet mounts to accept the standard monocular or binocular observation tubes (not included); both discussion and observation tubes independently rotatable 360° for the desired observation position	\$ 498.00
51	.3 266	Comparison tube with built-in prism system, half stops, adjustment screws to correct slight magnification differences of the two microscope objectives of equal image size and bayonet mount to accept standard microscope tube	\$ 640.00
		STAGES	
	12 028	Plain stage #19, with one pair of slide clips, 130 x 125mm	\$ 14.00
	51 321	Attachable mechanical stage SK1, traversing area 75 x 26mm	39.00
51	13 155	Attachable mechanical stage #21, with low set coaxial control knobs, traversing area 75 x 26mm	67.00
51	12 123	Mechanical stage #26, with scales and verniers, traversing area 76 x 50mm (3" x 2")	79.00
08	51 322	Mechanical stage #47, 135 x 135mm, with scales and verniers, low set operating knobs, traversing area $76 \times 50 \mathrm{mm}$	94.00
51	12 025	Mechanical stage #48, 160 x 138mm, with scales and verniers, low set operating knobs on one axis, traversing area $76 \times 50 \mathrm{mm}$ (3" x 2")	168.00
		BRIGHTFIELD CONDENSERS	
51	12 020	Single lens condenser #65, A.0.65, with iris diaphragm, filter holder with blue and ground glass filters, in cylindrical sleeve mount (only fitting the stand #051 230)	\$ 22.00
51	12 021	Two lens condenser #66, A. 1.20, with iris diaphragm, filter holder with blue and ground glass filters, in cylindrical sleeve mount (only fitting the stand #051 230)	36.00
0	51 342	Two lens condenser #70, A. 1.20, with iris diaphragm, filter holder with blue and ground glass filters, in dovetail mount	15.00
5	12 085	Swing-out condenser with aperture diaphragm and interchangeable aspherical top element N.A. 0.90, $\#601$	51.00
		PHASE CONTRAST CONDENSER	
5	13 140	Phase contrast condenser, PHACO #402a, with aperture diaphragm, swing-out upper element Achr. 0.90 and revolving disc containing three annular diaphragms, in dovetail mount	\$ 199.00
		DARKFIELD CONDENSERS	
5	13 109		\$ 91.00
5	13 111	Dry darkfield condenser #84, D 0.80, in centering mount	120.00
		ILLUMINATORS	
5	12 078	Plane and concave mirror12	\$ 12.00
5	14 096	Micro-dia lamp13, 15 watt, for direct connection to 110 volts A.C.	23.00
0	51 475	Low voltage lamp 6.3 volts, 0.8 amps, 5 watt with adjustable aspherical condenser, blue and ground glass filters and regulating transformer (mounted to the microscope by means of a bayonet lock)	56.00
5	14 095	•	85.00
		TRANSFORMERS	
22 5	500 124	Regulating transformer with voltmeter for 6 volts, 15 watts or 6 volts, 5 watts; for con-	
. •		nection to 110/120 volts, 60 cycles A.C.	\$ 67.00

	SPARE BULBS		
500 013	Bulb 15 watts, 110 volts A.C. (for #514 096)	\$	1.00
050 520	Bulb 6.3 volts, 5 watts (for #051 475)	•	0.40
500 012	Bulb 6 volts, 15 watts (for #514 095)		3,00
	CASE		
512 045			
312 045	Leatherette carrying case, upright	\$	34.00
	MISCELLANEOUS ACCESSORIES		
513 093	Diamond object marker	\$	42.00
513 030	Cover glass thickness gauge with graduation (1 interval = $1/100$ mm)	Ψ	31.00
513 211	Funnel stop for achromatic objective $63x/0.85$	\$	2.00
513 069	Funnel stop for achromatic objective $100x/1.30$ and fluorite $95x/1.32$		2.00
513 212	Funnel stop for fluorite $54x/0.95$		2.00
	NEW DRAWING ATTACHMENT		

NOTE:

microscopes

513 261

The microscope light source must be equipped with a regulating transformer. A 60 watt desk lamp is also required to illuminate the drawing area

Drawing attachment with $80/20~\mathrm{prism}$, adjustable lateral drawing tube, focusing device with built-in achromatic objective, eyepiece tube for the projection eyepiece (eyepiece not included), 45° mirror and bayonet mount for observation tubes; for SM and LABOLUX

RECOMMENDED EYEPIECES

519 051	Periplanatic widefield eyepiece, single GF 12.5x, field of view 18mm (drawing area 14 - 23cm)	\$ 49.00
519 137	Periplanatic widefield eyepiece, single GF 10x, field of view 18mm (drawing area 10 - 18cm)	39.00
519 180	Huygens eyepiece, single H 6.3x, field of view 18mm (drawing area 7 - 11cm)	16.00

Prices and Specifications are subject to change without notice, as we endeavor constantly to improve our products in keeping with the latest technical developments.

TERMS: NET 30 DAYS

Catalog #51.2-D.SM Printed in U.S.A.

\$ 301.00

Effective: January 1st, 1968



E. LEITZ INC. 468 PARK AVENUE SOUTH, NEW YORK, NEW YORK 10016

BRANCH OFFICES; SAN FRANCISCO, CALIF. AREA: 830 MITTEN ROAD • BURLINGAME, CALIFORNIA 94010 CHICAGO, ILL. AREA: 205 WEST TOUHY AVENUE • PARK RIDGE, ILLINOIS 60062

Catalog No. 51.2-D.SM 1/68/12M